

Services of the lab kitchen / scullery (“Medienküche”)

Autoclaving media:

- sterilisation is performed each Thursday
- weigh in media or agar composition without water, make a note with the necessary amount of water that has to be added, mention your name, and your room number
NOTE: Bring the bottles to the scullery until Wednesday afternoon!

General services:

- autoclaving of solutions, consumables like micro tubes or glassware and waste (standard: 121°C, 20 min)
- cleaning of lab dishes
NOTE: Bring used lab ware and liquid waste bottles directly to the scullery onto the designated red trays!
- refilling of 70% ethanol containers (e.g. below fume hood room. 2.28),
NOTE: Bring empty containers directly to the scullery do not put them back!

General lab equipment (glass bottles, pipettes, etc.):

equipment	storage location
autoclaved glass pipettes for cell culture	corridor locker next to room 2.28
bottles for liquid waste, signed with “Abfall” [only for S1 -liquid waste]	shelves on corridor next to room 2.29
Plastic bottles with a yellow cap [only for S2 -liquid waste]	room 2.22
S1 autoclave waste bin (w/o lid)	room 2.22
S2 autoclave waste bin (with lid and labelled with S2)	room 2.22
autoclaved Milli-Q water (desalted water); bottles with red caps	shelves on corridor next to room 2.29
sterile lab bottles (blue caps without any further mark)	metal cupboard on corridor next to room 2.05
sterile pipette tips for cell culture in boxes marked with red tape	corridor locker next to room 2.28
refill system for pipette tips– if you want to refill them yourself	corridor locker next to room 2.28