

PyRAT

Animal facility software








To familiarise yourself with PyRAT without risking errors in our own data, we recommend testing the demo version in advance. In the demo database, virtual animals can be selected, actions performed, work orders triggered, etc.




Under *Animals* > *Animals* all existing animals are displayed in an overview list.

1. View animals in the overview list






1. View animals in the overview list



2. Enter information for single animals


Ear punch numbers are entered here as *LabID*.

In the condition field, abnormalities such as "grooming" are entered..



3. Enter information for several animals at the same time


really is exactly the same action, i.e. the reason for killing or the condition are the same.



-5





4. View cages in the overview list



5. Create orders

(explained using the example of "Transport to the DRFZ")








5. Create orders

(explained using the example of "Transport to the DRFZ")

You can find your created requests and the processing status under Requests > Work requests.










Important: For new matings, imports of new strains or revitalisations, a final assessment must first be sent to tierarzt@drfz.de so that the new strain, associated mutations and further information can be entered into the PyRAT system. This also applies to orders placed with commercial breeders.

8. ET-Lab request creation




8. ET-Lab request creation






9. Virtual killing of pups

If animals are used for an experiment / organ removal before weaning




9. Virtual killing of pups

If animals are used for an experiment / organ removal before weaning





10. The animal pool

The animal pool includes animals that have been shared by scientists for use by others. These animals can be ordered and used by all scientists. If there is no call from the animal pool within the next 14 days, the request is considered a kill request with the reason "Surplus". Wild types are excluded - these remain in the animal pool indefinitely.




10. The animal pool





12. Change information after the death of the animals







You need help or have other concerns?

BeFor questions or other matters concerning PyRAT, please contact tierarzt@drfz.de

By the way, you can find a detailed PyRAT manual on the Intranet.



